



Cod	EventType	Sub	Inst/Con/Dis
b1	Basal Temp Rate/Start	Y	Disjunct
b9	Basal Temp Rate/End		End
p1	Basal Temp Percent/Start	Y	Disjunct
p9	Basal Temp Percent/End		End
bC	Basal Pattern Change	Y	Instant
p	Blood Pressure	Y	Instant
B	Bolus Delivery	Y	Instant
Bs	Bolus Square/End	Y	Instant
Bw	Bolus Wave/End	Y	Instant
C	Chloresterol	Y	Instant
X1	Exercise/Start	Y	Disjunct
X9	Exercise/End		End
F	Food Intake	Y	Instant
Gb	Glucose Blood	Y	Instant
Gs	Glucose Sensor	Y	Conjunct
H	Glycolated Hemoglobin	Y	Instant
I1	Insulin Active/Start	Y	Disjunct
I9	Insulin Active/End		End
IM	Insulin Manual	Y	Instant
K	KeyTone	Y	Instant
N	Note	Y	(Any)
PC	Pump Change	Y	Conjunct
IC	Pump Insulin Change	Y	Relative
IL	Pump Insulin Length	Y	Relative
Pp	Pump Prime	Y	Instant
Pr	Pump Rewind		Instant
P1	Pump Suspension/Start		Disjunct
P9	Pump Suspension/End		End
RC	Reservoir Change		Instant
S1	Sleep/Start		Instant
S9	Sleep/End		Instant
sC	Set Change		Instant
TC	Tubing Change		Instant
W	Weight	Y	Instant



Entity Type

- Reference
- Identifying
- Identifying/Major
- Event
- Event/Reading
- TransactionDetail
- History, Audit

Data Value