



EventType
 b Basal
 T1 Basal Temperature/Start
 T9 Basal Temperature/End
 P Blood Pressure
 B Bolus Reading
 C Chloolesterol
 X1 Exercise/Start
 X9 Exercise/End
 Gb Glucose Reading/Blood
 Gs Glucose Reading/Sensor
 HbA1C
 K KeyTone
 M Meal
 N Note
 Pr Pump Rewind
 Pp Pump Prime
 P1 Pump Suspended/Start
 P9 Pump Suspended/End
 W Weight

Open
 ReadingBasalTemperatureEnd ?
 ReadingBolusEnd ?
 Basal Pattern ?
 BasalTemp Pattern ?
 Pump ?
 Sensor ?

EventFlag
 E Error
 S Suspect
 D Doctor
 F Follow-up

GlucoseEventType
 M1 Pre Meal
 M9 Post Meal
 N Before Bed
 F Fast
 Ho Hypo
 H High
 Hr Hyper
 X1 Exercise/Before
 X9 xercise/After

BolusType
 N Normal
 S Square Wave
 D Dual Wave

MealType
 B Breakfast
 L Lunch
 D Dinner
 P Supper
 K Snack
 N Midnight