



**GlucoseEventType**

M0	Pre Meal
M1	1-Hour Post Meal
M2	2-Hour Post Meal
M3	3-Hour Post Meal
N0	Before Bed
N9	Overnight
F	Fast
X1	Exercise/Before
X9	Exercise/After

**BolusType**

N	Normal
S	Square Wave
D	Dual Wave

**EventType**

Cod	EventType	Sub	Inst/Con/Dis
b1	Basal Temp Rate/Start	Y	Disjunct
b9	Basal Temp Rate/End		End
p1	Basal Temp Percent/Start	Y	Disjunct
p9	Basal Temp Percent/End		End
cb	Basal Pattern Change	Y	Instant
P	Blood Pressure	Y	Instant
B	Bolus Delivery	Y	Instant
Bs	Bolus Square/End	Y	Instant
Bw	Bolus Wave/End	Y	Instant
C	Cholesterol	Y	Instant
X1	Exercise/Start	Y	Disjunct
X9	Exercise/End		End
M	Food Intake	Y	Instant
Gb	Glucose Reading/Blood	Y	Instant
Gs	Glucose Reading/Sensor	Y	Conjunct
H	Glycolated Hemoglobin	Y	Instant
I1	Insulin Active/Start	Y	Disjunct
I9	Insulin Active/End		End
IM	Insulin Manual	Y	Instant
K	KeyTone	Y	Instant
N	Note	Y	(Any)
Pc	Pump Change	Y	Conjunct
Pp	Pump Prime	Y	Instant
Pr	Pump Rewind		Instant
IP	Pump Insulin Length	Y	Relative
P1	Pump Suspended/Start		Disjunct
P9	Pump Suspended/End		End
CR	Reservoir Change		Instant
S1	Sleep/Start		Instant
S9	Sleep/End		Instant
CS	Set Change		Instant
CT	Tubing Change		Instant
W	Weight	Y	Instant

**MealType**

B	Breakfast
L	Lunch
D	Dinner
S	Second
P	Supper
K	Snack
T	Treatment
N	Midnight

**EventFlag**

D	Doctor
E	Error
F	Follow-up
S	Suspect

**Entity Type**

- Reference
- Major Identifying
- Identifying
- Transaction
- TransactionDetail
- History, Audit