

Damon DM Event V1.2



PumpHistory	
Name	
DateTime	
Model	
Type	
GlucoseTarget	
BasalRate	
CarbohydrateRatio	
CorrectionFactor	
InsulinEndDateTime	
Sensitivity	

Patient	
PatientId	
FirstName	
LastName	
BirthDate	
CarbTypeIsWeight	
IsImperialUnits	
EtCetera	

GlucoseEventType	
GlucoseTypeCode	
Description	

BolusType	
BolusTypeCode	
Description	

BasalPattern	
Name	
DateTime	
Rate	

CarbohydratePattern	
Name	
DateTime	
Amount	

EventFlag	
EventFlagCode	
Description	


EventType	
EventTypeCode	
Description	

Event	
DateTime	
EventTypeCode	
EventFlagCode	

- EventType
- b Basal Temp
 - P Blood Pressure
 - B Bolus Reading
 - BW Bolus Reading/Wave
 - C Chloresterol
 - X1 Exercise/Start
 - X9 Exercise/End
 - Gb Glucose Reading/Blood
 - Gs Glucose Reading/Sensor
 - g Glycolated Hemoglobin
 - K KeyTone
 - M Meal
 - N Note
 - Pr Pump Rewind
 - Pp Pump Prime
 - P1 Pump Suspended/Start
 - P9 Pump Suspended/End
 - W Weight

- EventFlag
- D Doctor
 - E Error
 - F Follow-up
 - S Suspect

Open
ReadingBasal is Temp
* Rate not %
* Need Start/End Rate ?



- GlucoseEventType
- M1 Pre Meal
 - M9 Post Meal
 - N Before Bed
 - F Fast
 - Ho Hypo
 - H High
 - Hr Hyper
 - X1 Exercise/Before
 - X9 xercise/After

- BolusType
- S Square Wave
 - D Dual Wave

- MealType
- B Breakfast
 - L Lunch
 - D Dinner
 - P Supper
 - K Snack
 - N Midnight